
Fig. 1

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The parallel plate capacitor is used to demonstrate fundamental electrostatic concepts and laws. The distance between the plates can be adjusted up to a maximum of 70mm, and it has a scale with millimeter graduations and a vernier which can be used to measure accurately down to 1/10mm.

Examples of experiments

- Experiments in a homogeneous electrical field (electrostatic inductance, separation of charge, transport of charges, force exerted on a charged body, image force)
- The relationship between charge, voltage and capacitance
- The effect of the distance between the plates and the dielectric material on capacitance
- Electrical field constant ϵ_0 , dielectric constant ϵ

References:

Description of experiments in the Main Catalogue of Physics Experiments, Volume 2 "Electricity - Electronics" (599 932)
"I-Measuring Amplifier D" (532 032)

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1 Safety Notes

- Use appropriate sources of voltage, which supply non-hazardous voltage levels. If no suitable power supply unit is available, limit the voltage output to a maximum of 2 mA using a safety fuse.

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2 Description, technical data

- ① main stand with two mounting surfaces base area: approx. 32cm X 15cm
(1.1) mounting surface for vertical arrangement of plates
(1.2) mounting surface for horizontal arrangement of plates
- ② fixed capacitor plate on insulator (2.1), with 4mm connection socket (2.2)
- ③ adjustable capacitor plate with 4mm connection socket (3.1), distance between plates can be adjusted to any point between 0 and 70 mm, approximate adjustment of the distance between the plates is made by sliding plate ③ with the stop screw ④ loosened, fine adjustments are made with using the adjusting knob ⑥, with the stop screw ④ tight
- ④ stop screw for the sliding arrangement for plate ③
- ⑤ scale with millimeter graduations and vernier for reading off the distance between the plates
measuring range: 70mm
measuring accuracy: 0.1mm

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- ⑥ adjusting screw for fine adjustment of the distance between the plates (actuate only with the stop screw ④ tight),
adjustment range: up to 20mm

Other technical data:

Diameter of the plates: 256mm

Thickness of the plates: 7mm

Capacitance with the plates 1mm apart: approx. 550pF

Deviations from the theoretical capacitance caused by marginal and leakage fields: approx. 40pF (average value for distances between plates of 1 and 40 mm.)

Mass: approx. 6.5kg

3 Use

Required in addition:

Voltage source, e.g.

stabilised power pack, 0...300 V DC 522 35
with measuring resistor, 10 M Ω (for limiting current) .. 536 21

or

high voltage power pack, 10 kV 522 37

or

glass tube 541 01

plastic rod 541 04

Material for rubbing (wool cloth, dry newspaper)

Static voltage measuring device ($R_i \geq 10^9 \Omega$), e.g.

static voltmeter 1.5kV 540 38

or

electrometer amplifier 532 14

with STE capacitors 1nF, 10nF, 100nF 578 25/10/31

voltage measuring device,

measuring range 1 V - 10 V DC 531 91

Device for measuring charge e.g.

I-measuring amplifier D 532 00

or

electrometer amplifier 532 14

with plug-in capacitors 1nF, 10nF, 100 nF 578 25/10/31

voltmeter, measuring range

1 V - 10 V DC for amplifier voltage, e.g. 531 91

High voltage cable (for voltages ≥ 1 kV 501 05

Earth connection for the experimental apparatus, e.g.

Connecting rod 532 16

Experimental cable, 1m, e.g. 501 42

Recommended if the relative humidity of the air is high:

Warm air fan

When carrying out experiments on the laws of plate capacitors, as well as when determining ϵ_0 and ϵ , work with the smallest possible distance between the plates, (interference caused by the marginal field increases with increasing distance between the plates). If necessary, allow for additional capacitances caused by the marginal and leakage fields as well as cables and measuring instruments.

First, with the stop screw ④ loosened, make approximate adjustment to the distance between the plates by sliding. Make the fine adjustment with stop screw ④ tightened.

Make sure the apparatus is properly earthed, (earth as a fixed reference potential). Earth the apparatus and, if necessary, the earth socket.

If the relative humidity of the air is high, it is advisable to blow warm air between the plates, to prevent the capacitor from discharging rapidly.

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Fig. 2
Measurement of the capacitor discharge as a function of the charge voltage and of the distance between the plates measuring the charge using an I-measuring amplifier D (532 00) with a voltage measuring device connected

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Fig. 3
Measurement of the voltage between the plates of a charged plate capacitor as a function of the distance between the plates
Quasi static measurement of voltage via a capacitive voltage separator at the input to the electrometer amplifier (532 14), with a voltage measuring device connected

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Fig. 4
Influence, separation of charge
Measurement of the charge of each influence plate (546 13) with the electrometer amplifier switched as a charge measuring device (532 14), with a voltage measuring device connected

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Fig. 5
Measurement of the force which is encountered by a charged experimental body in the (nearly) homogenous electrical field of the plate capacitor
Measurement of the force, using a DMS* force distance recorder 0.2 N (314 26) and a Newton meter (314 25) - Fig. 5.1 -
or with current balance (516 32) - Fig. 5.2

*) DMS expansion measuring strip